

## New Products

### Actimove® ManusForte – Functional Wrist Brace

Made from top quality, latex-free materials Actimove® Manus Forte wrist brace is used for the treatment of acute or chronic pain and irritation in the wrist joint, tendinous synovitis, arthrosis and postoperative or posttraumatic irritation.


The wrist brace ensures strong support and comfortable stabilisation for effective pain relief and advanced recovery. Actimove® Manus Forte leaves the fingers completely unrestrained, enabling the patient to take part in everyday activities.

<ul style="list-style-type: none"> <li>• Wide inelastic straps for extra strong stabilisation</li> <li>• Anatomically shaped orthopaedic splint prevents flexion and extension of the wrist</li> <li>• Two plastic stays on the backside support stabilization of the wrist</li> <li>• Safe hook and loop fasteners for easy handling and strong fixation</li> <li>• Highly breathable, skin-friendly materials provide high wearing comfort</li> <li>• Self-explanatory, tubular construction for easy application with one hand</li> </ul>	 <p>Enlarge photo</p>
--	--

### Actimove® FootLift - Stabilising Leaf Spring AFO

Actimove® FootLift is a specially designed ankle foot orthosis made of polypropylene and fabricated in a predorsiflexed position. It enables physiological walking and reduces the risk of stumbling. Actimove® Footlift is indicated for peroneal paralysis as well as drop foot secondary to neurological pathologies caused by stroke, cerebral palsy or trauma.



 <p>Enlarge photo</p>	<ul style="list-style-type: none"> <li>• Anatomically shaped foot plate</li> <li>• Adjustable calf band fits a broad range of patients</li> <li>• Heat mouldable for individual fit</li> <li>• Non-slip, soft calf padding provides good cushioning</li> <li>• Easy to apply and remove</li> <li>• Slim design for a good fit, even in small shoes</li> </ul>
--	---

### Actimove® EpiSport –Anatomically Shaped Elbow Brace

Actimove® EpiSport is a specially designed elbow brace made from elastic, flexible, durable and latex free material. It is indicated for any kind of acute and chronic Epicondylitis e.g. associated with playing golf or tennis.

 <p>Enlarge photo</p>	<ul style="list-style-type: none"> <li>• Targeted dual compression relieves pressure of affected muscles, for effective pain relief</li> <li>• Anatomically correct design for comfortable fit</li> <li>• Buckle closure and Velcro® Strap for easy handling</li> <li>• Elastic strap prevents blood constriction</li> <li>• Non-slip and soft cushioning inner liner ensures optimal fit</li> </ul>
--	--

### Actimove® StepLite – Anatomically shaped heel cushions

Actimove® StepLite heel cushions are made of medical grade, odourless silicone. They are indicated for the treatment of plantar fasciitis, calcaneal spur and pain in the Achilles tendon. They can also be used for the compensation of a difference in leg lengths of 1 cm.

The heel cushions with an embedded BlueSpot for the most sensitive area of the heel efficiently reduce pressure and pain. Actimove® StepLite has been specially designed to give an exceptional feeling of stabilisation and comfort ensuring improved patient compliance and a smooth healing process.

 <p>Enlarge photo</p>	<ul style="list-style-type: none"> <li>• Effective pain relief due to embedded BlueSpot</li> <li>• Targeted shock absorption</li> <li>• Sleek anatomical shape fits neatly into most shoes</li> <li>• Good bonding inside the shoe to ensure a secure fit</li> <li>• Non-slip surface for secure fit</li> <li>• Medical grade, odourless silicone reduces the risk of allergic reactions</li> </ul>
---	---

### Actimove® TaloCast-Air - Functional Ankle Brace

Actimove® TaloCast-Air is a functional ankle brace indicated for acute capsular ligament injuries, ankle strains and chronic ankle instability. It can also be used for preventative bracing and post cast support.

With its inflatable cushioning air cells the ankle brace stabilizes the joint comfortably giving patients a feeling of security. The sleek anatomic design allows patients to move freely after injury for a quick recovery.

<ul style="list-style-type: none"> <li>• Individually adjustable air cells create a massage effect helping to reduce swelling</li> <li>• Semi-rigid, lightweight shells give the right balance of stability and mobility</li> </ul>	
---	--

- High-quality durable heel pad and 360° surround straps provide a strong and tight fit
- Heel pad can be easily adjusted with hook and loop closure
- The heel pad`s printed foot symbol ensures correct application
- Comfortable fit in everyday footwear due to sleek and slim design



Enlarge photo

### Actimove® TaloCast-AirGel - Functional Ankle Brace

Actimove® TaloCast-AirGel is a functional ankle brace indicated for acute capsular ligament injuries, ankle strains, chronic ankle instability as well as swelling and edema at the ankle. It can also be used for preventative bracing and post cast support.

With its inflatable cushioning air cells the ankle brace stabilizes the joint comfortably giving patients a feeling of security. The sleek anatomic design allows patients to move freely after injury for a quick recovery. Two integrated gel pads allow the application of cryotherapy for efficient pain reduction.



Enlarge photo


- Integrated cooling gel pads effectively help to reduce swelling, edema and pain, promoting the healing process
- Individually adjustable air cells create a massage effect further helping to reduce swelling
- Gel pads and aircells comfortably cushion and protect the ankle joint
- Semi-rigid, lightweight shells with gel pads provide increased stabilisation
- High-quality durable heel pad and 360° surround straps provide a strong and tight fit

	<ul style="list-style-type: none"> <li>• Heel pad can be easily adjusted with hook and loop closure</li> <li>• The heel pad`s printed foot symbol ensures correct application</li> <li>• Comfortable fit in everyday footwear due to sleek and slim design</li> </ul>
--	---

### Delta-Lite® Plus - The Trusted Rigid Solution.

Delta-Lite® Plus is our new fiberglass cast tape that combines a fiberglass substrate and advanced water activated polyurethane resin system which provides tack free properties. The tack free properties allow ease of application and molding using standard examination gloves.

Unique Features:

<ul style="list-style-type: none"> <li>• New resin formulation -use exam glove of your choice</li> <li>• Rigid support</li> <li>• High lamination</li> <li>• Smooth surface</li> <li>• Weight bearing in approximately 20 minutes</li> <li>• 3-5 minute set time</li> </ul>	 <p data-bbox="732 1503 903 1532">Enlarge photo</p>
---	---


### Delta-Cast®Soft - The Flexible Solution.

Delta-Cast® Soft is a semi-rigid non-fiberglass cast tape, ideal for customized support. Delta-Cast® Soft provides variable support for fractures and soft tissue injuries while allowing effective stabilization for the injured area. The use of a semi-rigid cast material allows for greater range of movement, earlier mobility and reduces some of the problems of rigid immobilization, such as muscle atrophy.

Indications of use:

- Semi-rigid support
- Effective stabilization/Functional immobilization

Unique features:

<ul style="list-style-type: none"> <li>• Versatile application</li> <li>• Good conformability</li> <li>• Rapid removal with scissors or by unwrapping</li> <li>• Patient comfort</li> </ul> <p>Smooth finish Soft edges</p> <ul style="list-style-type: none"> <li>• Excellent for circular and removable casts</li> </ul>	 <p>Enlarge photo</p>
--	--

## Cutimed® Sorbact® gel

### The antimicrobial dressing with hydrogel

The new ready-to-use hydrogel impregnated dressings reduce bacterial load by combining the unique antimicrobial action of Cutimed Sorbact with moist wound management.

The hydrogel component is free of preservatives and promotes autolytic debridement in sloughy or partly necrotic wounds.

<ul style="list-style-type: none"> <li>• No risks or side effects compared to the use of antiseptics or local antibiotics</li> <li>• No development of bacterial or</li> </ul>	
--	--

- fungal resistance
- No cytotoxicity
- Supports the natural wound healing process



[Enlarge photo](#)

[Click to learn more about the Cutimed Sorbact range.](#)

## **Delta-Dry® Water Resistant Cast Padding**

The solution that fits your life.

Delta-Dry®'s unique patented technology outperforms other water resistant cast padding, making it the best choice for clinicians and a convenient solution for patients over traditional cast padding.



[Enlarge photo](#)

- Allows patients to shower, swim and bathe like normal
- No special tools are needed for application or removal
- Conformable enabling patient to having increased comfort and protection
- Easy to apply
- Better long term water resistance
- Offers more under cast protection
- Eliminates the need for stockinette

--	--

## **JOBST® Elvarex® Soft**

Ready-To-Wear Flat Knit Medical Compression Stockings

Keeping the high standards of JOBST® Elvarex® in mind, our engineers have developed a new softer textile for the management of mild to moderate edema. The JOBST Elvarex Soft Lower Extremity Stockings have been designed to apply a moderate compression while making it easy to put on and take off.



[Enlarge photo](#)