

Leukomed® Sorbact®

Cover dressing coated with DACC Technology



Leukomed® Sorbact®

Cover dressing coated with DACC Technology


This innovative cover dressing not only helps to protect a wound from external contamination, but due to its contact layer coated with DACC Technology, helps reduce colonization of harmful microbes.

Adhesive film

- **Breathable and shower-proof**
Provides high patient comfort
- **Bacterial barrier**
Protects against external contamination
- **Reliable adhesion**
Dressing stays securely in place

Absorbent wound pad

- **Absorbs low - moderate fluid levels**
Reduces the risk of maceration





Wound contact layer coated with DACC Technology

- **Safely¹ and effectively binds bacteria and fungi without a chemically active agent**
Helps to minimize the risk of infection²

Indications:

Wounds with low to moderate exudate such as:

- Surgical incisions
- Post-operative dehisced wounds
- Traumatic wounds
- Lacerations, cuts, abrasions
- Minor burns



Order information

Size	Pad Size	REF No.	Dressings/Box	Case Contents	Suggested HCPCS*
2 in x 3 in	1.1 in x 1.5 in	7619900	20	56	A6251
3 in x 4 in	1.6 in x 2.6 in	7619901	20	12	A6203
3 in x 6 in	1.6 in x 4.3 in	7619902	20	8	A6203
4 in x 8 in	2 in x 6 in	7619903	20	12	A6203
4 in x 10 in	2 in x 8 in	7619904	20	10	A6204
4 in x 12 in	2 in x 10 in	7619905	20	8	A6204
4 in x 14 in	2 in x 12 in	7619906	20	8	A6204

1. Haycocks S. et al. Use of a DACC-coated antimicrobial dressing in people with diabetes and a history of foot ulceration. Wounds UK (2011) 7: 11, 108-114.

2. Grothier L. Managing wound infection decisions via a color-coded patient pathway for improved patient outcomes and reduced healthcare costs. SAWC Poster Presentation 2012

Please refer to the product label and / or package insert for full instructions on the safe use of these products.

* This document includes codes that might be used to bill for BSN medical products and related services. It is the provider's sole responsibility to determine and submit appropriate codes, charges, and modifiers for services rendered. BSN medical Inc. cannot guarantee coverage or reimbursement.

BSN medical Inc. 5825 Carnegie Blvd. Charlotte, NC 28209-4633
Tel. (+1) 704 554 9933 Fax (+1) 704 358 4558 www.bsnmedical.com
To order toll-free: BSN medical (+1) 800 552 1157

